



Video Laryngoscope & Nishikawa Blade

Perform Routine and Difficult Tracheal Intubation with all Confidence.

Nishikawa Blade

S-shaped Spatula


Designed to manage the enlarged tongue base.

Bifid Tip

Fits better into vallecula space & lift epiglottis more effectively.

S-shaped spatula and bifid tip

Designed to manage the bulk of a large tongue.



(From Nishikawa K, Yamada K, Sakamoto A. A new curved laryngoscope blade for routine and difficult tracheal intubation. Anesth Analg. 2008;107:1248-1252 with permission)

E.AN-II PLUS

3" High Resolution Full Angle

16:9 Screen

Rechargeable Battery

Water Proof (IP64)


Nishikawa Blade


Camera & Video

Anti Fog

Metal Handle

High Transparency Blade






Nishikawa Laryngoscope Blade

This blade has an S-shaped Spatula and Bifid Tip, designed to manage the bulk of a large tongue.

From Nishikawa H, Yamada K, Sakamoto A. A new curved laryngoscope blade for routine and difficult tracheal intubation. Anesth Analg. 208; 107:1248-1252 with permission.

Quote from ATLAS OF AIRWAY MANAGEMENT
TECHNIQUES AND TOOLS - P33



E.AN-IIL - E.AN-IIS

MODEL	SIZE	PACKAGE	BEST PATIENT WEIGHT
E.AN-IIL	Standard Adult	1	40-120 kg
E.AN-IIS	Standard Pediatric	1	1.5-50 kg




E.AN-IILP - E.AN-IISP

MODEL	SIZE	PACKAGE	BEST PATIENT WEIGHT
E.AN-IILP	Premium Adult	1	40-120 kg
E.AN-IISP	Premium Pediatric	1	1.5-50 kg



NISHIKAWA BLADES

MODEL	SIZE	PACKAGE	BEST PATIENT WEIGHT
ZH-LS	114mm/16mm	200/case-50/box	50-120 kg
ZH-MS	106mm/15mm	200/case-50/box	40-100 kg
ZH-SS	83mm/11mm	200/case-50/box	4-50 kg



REGULAR BLADES

MODEL	SIZE	PACKAGE	BEST PATIENT WEIGHT
ZH-L	114mm/16mm	200/case-50/box	50-120 kg
ZH-M	106mm/15mm	200/case-50/box	40-100 kg
ZH-S	83mm/14mm	200/case-50/box	4-50 kg
ZH-N	79mm/9mm	200/case-50/box	1.5-5 kg



CHANNEL BLADES

MODEL	SIZE	PACKAGE	MAX. TRACHEAL TUBE SIZE BEST PATIENT WEIGHT
ZH-LG	114mm/20mm	100/case-25/box	Ø 8.5 50-120 kg
ZH-MG	106mm/19mm	100/case-25/box	Ø 7.5 40-100 kg